



DONJOY[®]
VERSAROM[™] ***HIP***

PRECISE, INDIVIDUALIZED
PROTECTION AND SUPPORT.



PRECISE, INDIVIDUALIZED PROTECTION AND SUPPORT.

The low-profile, DonJoy® VersaROM™ Hip helps deliver effective, functional ROM restriction in 15° increments. Easy-to-use abduction/adduction push buttons and malleable hip and thigh shells help provide an individualized fit for increased patient comfort and compliance. The adjustable compression belt with easy-grip pull tabs anchors the hip joint, and the padded thigh cuff increases optimal alignment. With a universal fit, one brace fits right or left hip up to a 137cm (54") waist circumference.

INDICATIONS: DonJoy® VersaROM™ Hip is indicated for use with post-op patients who have undergone a surgical procedure and need an orthosis until the component is integrated and soft tissue heals.

ADDITIONAL INDICATIONS INCLUDE:

- Hip instability
- Joint revisions or reconstruction
- Stable hip fractures
- ROM control
- After reposition of luxated hip joints
- Muscular imbalance
- Arthroscopic hip procedures
- Labral repair, with or without gluteus medius repair
- Post-op proximal hamstring repair
- Patients at risk to dislocate



**INTEGRATED BELT COMPRESSION
WITH CINCHING STRAPS:**

Pull tabs anchor the hip joint
and thigh cuff to promote
optimal alignment

**MALLEABLE HIP AND
THIGH FRAMES:**

Conforms to patient's hips
for optimal, individualized fit

PRECISE ROM ADJUSTMENTS:

Range of Motion allows
controlled hip flexion
and extension adjustments
in 15° increments

**TELESCOPING HEIGHT
ADJUSTMENT:**

Helps to accommodate a
wide range of heights and
body shapes

**LATERALLY PLACED
ABDUCTION PADDLE WITH
THREE-POINT PRESSURE:**

Helps to maintain the head
of the femur in the acetabulum

CUSTOMIZED FIT:

Two optional push-button hinges
adjust abduction and adduction for
an individualized fit

**LIGHTWEIGHT, SLIM-
PROFILE DESIGN:**

Combined with the
softgoods, helps to maximize
patient comfort and
compliance

ORDERING INFORMATION

Part number	Description
11-3335	DonJoy® VersaROM™ Hip



DJO | 1a Guildford Business Park | Guildford | Surrey | GU2 8XG | UNITED KINGDOM

www.DJGlobal.eu

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.