



DONJOY[®]
VERSAROM[™] HIP

VersaROM™ Hip



Integrated belt compression with cinching straps:
Pull tabs anchor the hip joint and thigh cuff to promote optimal alignment

Malleable hip and thigh frames:
Conforms to patient's hips for optimal, individualized fit

Precise ROM Adjustments:

Range of Motion allows controlled hip flexion and extension adjustments in 15° increments

Telescoping Height Adjustment:
Helps to accommodate a wide range of heights and body shapes

Laterally placed abduction paddle with three-point pressure:
Helps to maintain the head of the femur in the acetabulum

Customized Fit: Two optional push-button hinges adjust abduction and adduction for an individualized fit

Lightweight, slim-profile design: Combined with the softgoods, helps to maximize patient comfort and compliance

Not made with natural rubber latex

**For more information, please contact the
DJO Bracing & Supports Representative listed below:**



DJO, LLC
T 800.336.6569 | D 760.727.1280 | F 800.936.6569
2900 Lake Vista Drive | Lewisville, TX 75067 | U.S.A.
DJOglobal.com

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.